CERVOS
MEDICAL

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ENGLISH / US

Sterile: Contents sterile unless package is opened, damaged or expired.

DESCRIPTION: The Cervos Platelet Separator is composed of up to 3 process disposables, a threaded rod and up to 2 luer lock connectors.

- A. INTENDED PURPOSE:** The Cervos Platelet Separator is designed for the preparation of autologous platelet rich plasma (PRP) from a small sample of peripheral blood at the patient's point of care.
- B. INTENDED USER:** The Cervos Platelet Separator is intended to be used by a licensed physician familiar with the possible side effects, typical findings, limitations, indications, and contraindications of autologous platelet rich plasma.
- C. INTENDED PATIENT POPULATION:** The Cervos Platelet Separator is intended to be used on any patients that are suitable for autologous platelet rich plasma as determined by the physician.
- D. INDICATION(S):** The Cervos Platelet Separator is indicated for preparing autologous platelet rich plasma.
- E. CONTRAINDICATION(S):** The use of the Cervos Platelet Separator may be contraindicated when there is: Clinical or laboratory evidence of septicemia for patients who have taken aspirin, or other medications that alter platelet function, within 3 days prior to surgery, or patients with disorders associated with platelet dysfunction.
- F. INTENDED LIFETIME:** The Cervos Platelet Separator is a sterile, single-use device intended for
 - Transient use (<60 minutes)
- G. INTENDED CLINICAL BENEFIT:** The Cervos Platelet Separator is intended to prepare autologous platelet rich plasma (PRP) from a small sample of peripheral blood at the patient's point of care to be used as determined by the physician.

WARNINGS/CAUTIONS: For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient. There are no limitations of the device as determined by the manufacturer since there are no clinical claims being made for the device. Do not use if package is open or damaged and if the expiry date has been exceeded. Do not continue to use if the device gets damaged during procedure. Follow manufacturer instructions when using centrifuge. Use a programmable swinging bucket centrifuge capable of creating a sustained G force of 2300 for 2 minutes and 5 minutes respectively. Use a centrifuge with buckets that can fit a disposable with an outer diameter of 1.64" and a height of 4.62". Aseptic technique should be used during all processing steps including when aspirating and dispensing fluids. Use prepared PRP within 4 hours. Any unused PRP is to be discarded in accordance with local biohazard requirements. The PRP prepared by this device is not indicated for delivery to the patient's circulatory system. Possible adverse effects include: Damage to blood vessels, hematoma, delayed wound healing, temporary or permanent nerve damage that may result in pain or numbness and early or late post-operative infection are associated with surgical procedure. It is recommended that in order to counteract the citrate anticoagulant, 10% Calcium Chloride Injection, USP solution be added to the PRP in a ratio of PRP to Calcium Chloride solution of 20:1 prior to mixing the PRP with the bone graft material. After use, this product may be a potential biohazard. Dispose in accordance with applicable laws and regulations.

Note: These instructions are NOT meant to define or suggest any medical or surgical technique. The individual practitioner is responsible for the proper procedure and techniques to be used with this device. Use only a luer lock syringe with this device.

H. DIRECTIONS FOR USE

1. Swab luer port on the top of the separator with alcohol prior to each use.
2. Ensure that only blood from only one patient is processed. Attach a

syringe to the separator canister. Pull a vacuum by drawing to the end of the syringe and unscrew it while holding the plunger. Repeat this step for all separator canisters.

3. Attach syringe containing blood to luer port on the top of a grey separator canister and load approximately 10mL of blood, remove 5mL of blood to purge the air. Then slowly add anti-coagulated blood to either one or more grey separator canisters. Remove thumb from the syringe plunger to allow the pressure to equalize in the disposable prior to disengaging the syringe.
4. If using an even number of gray disposables, ensure that **equal amounts** of anti-coagulated blood are added to each disposable and only grey canisters are utilized. Ensure swabable connector is tight after filling disposable. Place loaded separator into the centrifuge.
5. If an odd number of separators are being used, use a counterbalance filled with an equal amount of water. Place the counterbalance in the centrifuge directly opposite from the separator.
6. Follow centrifuge manufacturer's instructions on setting the G force such that the separator is exposed to 2300 G's for 2 minutes with a 90 second minimum brake period.
7. Remove each disposable separator from the centrifuge and attach a male to female luer lock connector by pushing it into the connector and turning it to the right. Using a syringe and threaded rod, extract the plasma and platelets leaving behind the red cells.
8. Attach syringe containing ACD-A to the remaining separators, add 2mL of ACD-A and then connect the syringe containing plasma & platelets to the luer port on the top of a maroon separator canister and pull back on the plunger slightly. Then slowly add anti-coagulated platelets and plasma into one or more maroon separator canisters. Remove thumb from the syringe plunger to allow the pressure to equalize in the disposable prior to disengaging the syringe.
9. If using multiple canisters, ensure that **equal amounts** of anti-coagulated platelets and plasma are added to all canisters and only maroon canisters are utilized for the second centrifugation step. Ensure swabable connector is tight after filling disposable. Place loaded separator into the centrifuge.
10. If an odd number of separators are being used, use a counterbalance filled with an equal amount of water.
11. Follow centrifuge manufacturer's instructions on setting the G force such that the separator is exposed to 2300 G's for 7 and a half minutes with a 90 second minimum brake period. If using a centrifuge capable of creating a sustained G force of 2800, set the G force such that the separator is exposed to 2800 G's for 6 minutes with a 90 second minimum brake period.
12. Remove each separator from the centrifuge. Using a syringe and threaded rod extract 1mL of platelet poor plasma (PPP) into a syringe from one separator and set aside. Using another syringe and the threaded rod, extract the plasma leaving behind approximately 2mL-10mL of PRP in the separator.
13. Using a syringe, flush back and forth the remaining PRP to ensure thorough mixing and then pull a vacuum by drawing to the end of the syringe and unscrew it while holding the plunger. Using the syringe containing the 1mL of platelet poor plasma, flush back and forth each empty separator to remove any platelets that may have adhered to the surface of the separator. The two syringes of PRP can be combined or administered separately.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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Ranfac® Corp.