# Centrifuge CM-8S

User manual

## Contenuti

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### 1.1 Specifiche tecniche

Rotor rotation speed, RPM	from 100 to 4500
Relative centrifugal force, RCF (G)	from 1 to 3735
Accuracy of rotation speed maintenance, %	±2
Timer range	1 min. - 12 h
Rotor speed step (RPM)	1
Centrifugal force step, RCF (G)	1
Timer step	1 second
Number of centrifuge brake levels	10
Noise level (1 m distance, dBA)	no more than 55
Electromagnetic compatibility	according to EN 61326-1: 2013
Ambient conditions:	
Ambient temperature, °C	from +10 to +40
Relative humidity, %	80
Do not use in humid premises (above 80% humidity)	
Use up to a height of 2000 m above sea level	
Total test tube imbalance limit, g	Refer to pages 16, 17
Mains voltage, V	110-240 V
Mains frequency, Hz	50-60 Hz
Overvoltage category	II
Power, W	320
Size of the device (length x width x height), mm	431x410x247
Approximate weight, kg (without rotor and cord)	14.7
Safety area	300 mm
Max. permitted density of material to be centrifuged	1.2 g/ml
Pollution degree	2
Number of programs, created by user	10

### 1.2 Content of package

<b>Name</b>	<b>pcs</b>
Centrifuge	1
Power cord	1
Rotor	1
Rotor key / emergency lid opening key	1
Rotor nut and spacer	1
User manual	1
Packaging	1

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## 7 Rotor installation and usage

**Installation:** Place the rotor on the axis. On top of the rotor place a special spacer. Note that the pin of the spacer must get into the groove of the axis. Place the rotor nut on the axis and tighten it with the included rotor key. Before each start check that the rotor is firmly tightened. Disassembly is made in reverse sequence.

**Usage:** Always load the rotor symmetrically. Each test tube should be counterbalanced by another test tube. By minimizing the difference of weight between the filled test tubes, you reduce the wear of mechanical parts of the centrifuge. If the overall imbalance of test tubes in the rotor exceeds allowed value, centrifuge will start emergency braking and the display will show Erdb – error code for imbalance.


**Do not use damaged rotors!**


**ATTENTION!** Samples may contain pathological material, including pathogens of serious diseases. **Always check test tubes for damage before centrifugation, to prevent centrifuge contamination.**

**Attention!** To avoid imbalance of the rotor, all adapters should be of equal weight when loaded (according to the marking on the adapter, see the picture below). If this requirement is not observed, a strong imbalance of the rotor can occur. The weight is indicated on the adapter, see the picture below.

## 8 Applied rotors

Rotor	Description
	<p><b>ROTOR 8M</b> for 12 adapters.            Maximal volume of applied test tubes: 12 ml.            Maximal size of applied test tubes (LxD): 115 x 17.5 mm.            Top speed: 4500 rpm.            Total test tube imbalance limit: 5 g.            Weight: 730 g.  <b>Autoclave sterilization at a temperature up to 121 °C (20 minutes) is allowed.</b></p>
	<p><b>ROTOR 8M.01</b> for 4 adapters.            Maximal volume of applied test tubes: 50 ml.            Maximal size of applied test tubes (LxD): 130 x 30 mm.            It is possible to order inserts in 50 ml adapters (diameter 18 and 22 mm) for the test tubes (DxL) 17.5 x 115 mm/20 x 115 mm.            Top speed: 4500 rpm.            Total test tube imbalance limit: 3 g.            Weight: 940 g.  <b>Autoclave sterilization at a temperature up to 121 °C (20 minutes) is allowed.</b></p>
	<p><b>ROTOR 8M.02</b> for 24 adapters.            Maximal volume of applied test tubes: 12 ml.            Maximal size of applied test tubes:            Inward row: (D x L): 17,5 x 115 mm            Outward row: (D x L): 17,5 x 140 mm            Top speed: 4500 rpm.            Total test tube imbalance limit: 5 g.            Weight: 1065 g.  <b>Outer row is loaded first (test tubes L=140 mm, volume 15 ml)</b>  <b>Autoclave sterilization at a temperature up to 121 °C (20 minutes) is allowed.</b></p>
	<p><b>ROTOR 8M.03</b> for 4 adapters.            Maximal size of applied test tubes (LxD): 125 x 44 mm.            Top speed: 4500 rpm.            Total test tube imbalance limit: 5 g.            Weight: 1120 g.  <b>Autoclave sterilization at a temperature up to 121 °C (20 minutes) is allowed.</b></p>
	<p><b>ROTOR 8M.05</b> for 12 adapters with a working angle 10°.            Maximal volume of applied test tubes: 15 ml.            Maximal size of applied test tubes (LxD): 140 x 17.5 mm.            Top speed: 4500 rpm.            Total test tube imbalance limit: 5 g.            Weight: 800 g.  <b>Autoclave sterilization at a temperature up to 121 °C (20 minutes) is allowed.</b></p>
	<p><b>ROTOR 8M.06</b> for 6 adapters.            Maximal volume of applied test tubes: 50 ml.            Maximal size of applied test tubes (LxD): 120 x 30 mm.            It is possible to order inserts in 50 ml adapters (diameter 18 and 22 mm) for the test tubes (DxL) 17.5 x 115 mm/20 x 115 mm.            Top speed: 4500 rpm.            Total test tube imbalance limit: 5 g.            Weight: 1015 g.  <b>Autoclave sterilization at a temperature up to 121 °C (20 minutes) is allowed.</b></p>

Rotor	Description
	<p><b>ROTOR 8M.07</b>            Universal rotor for 4 anodized aluminum adapters (flat bottom), also suitable for platelet-rich plasma (PRP) preparation.            Maximal volume of applied test tubes: 50 ml.            Maximal size of applied test-tubes (D x L): 30 x 130 mm.            Top speed: 4500 rpm.            Total test tube imbalance limit: 5 g  <b>Autoclave sterilization at a temperature up to 121 °C (20 minutes) is allowed.</b></p>
	<p><b>ROTOR 8M.08</b>            Universal rotor for 6 anodized aluminum adapters (flat bottom), also suitable for platelet-rich plasma (PRP) preparation.            Maximal volume of applied test tubes: 50 ml.            Maximal size of applied test-tubes (D x L): 30 x 120 mm.            Top speed: 4500 rpm.            Total test tube imbalance limit: 5 g.  <b>Autoclave sterilization at a temperature up to 121 °C (20 minutes) is allowed.</b></p>
	<p><b>ROTOR 8M.09</b>            Universal rotor for 2 anodized aluminum adapters (cone bottom) and 4 stainless steel adapters.            Maximal size of applied test-tubes (D x L):            30 x 130 mm - 2 pcs.            17.5 x 120 mm - 4 pcs.            It is possible to order inserts in 50 ml adapters (diameter 18 and 22 mm) for the test tubes (DxL) 17.5 x 115 mm/20 x 115 mm.            Top speed: 4500 rpm.            Total test tube imbalance limit: 3g.  <b>Autoclave sterilization at a temperature up to 121 °C (20 minutes) is allowed.</b></p>
	<p><b>ROTOR 8M.10</b>            Universal rotor for 2 anodized aluminum adapters (flat bottom) and 4 stainless steel adapters, also suitable for platelet-rich plasma (PRP) preparation.            Maximal size of applied test-tubes (D x L):            37 x 125 mm - 2 pcs.            17.5 x 120 mm - 4 pcs.            Top speed: 4500 rpm.            Total test tube imbalance limit: 3 g.  <b>Autoclave sterilization at a temperature up to 121 °C (20 minutes) is allowed.</b></p>
	<p><b>ROTOR 8M.18</b>            Universal rotor for 2 anodized aluminum (cone bottom) and 2 stainless steel adapters.            Maximal size of applied test-tubes (DxL):            30 x 130 mm – 2 pcs.            17,5 x 131 mm – 2 pcs.            Top speed: 4500 rpm.            Total test tube imbalance limit: 3 g.  <b>Autoclave sterilization at a temperature up to 121 °C (20 minutes) is allowed.</b></p>

## 9 Troubleshooting (Error code table)

Error code	Sound signal	Cause	Correction
	-	No mains power connection. Power failure. No display on screen	Check connection of the power cable
	+	Internal problem	Contact local service office
	+	Total imbalance of test tubes is more than permitted	It is necessary to insert test tubes symmetrically (on the opposite sides). Fill in the test tubes up to the same weight. Each test tube should be counterbalanced by another test tube. Refer to pages 16, 17.
	+	Centrifuge is overheated	Check the air vents of device and give a little time and centrifuge will cool down automatically.
	-	Residual rotation of the rotor due mains power failure	Wait until the rotor is motionless; make sure there is not power failure.
	+	Centrifuge lid is left open	Close and press on the lid till you hear a click.
	+	One of control buttons ("OPEN"; "START/STOP") is stuck	Check the buttons for stuck. If the problem remains, contact local service office.

Microprocessor control system self-tests and reports detected faults and errors to the user.

## 10 Maintenance of centrifuge elements

Laboratory staff must pay special attention to the fact that the centrifuge parts of key importance are not damaged due to safety reasons.

Routine check	Frequency of checks	Not acceptable	Correction
Housing, plug, cord, control system, adapters (external inspection)	After every 160 hours of operation*	Cracks, compression marks, broken coating parts	Replace damaged parts
Condition of rubber parts: engine shock absorbers, instrument legs, sealing gum	Every two years or after every 3600 hours of operation	Cracks, increased rubber hardness	Replace damaged parts
Bearing condition	Every two years or after every 3600 hours of operation	Free motion is more than 0,05 mm, rotation of a rotor with jamming	Replace engine
Engine condition	Every two years or after every 3600 hours of operation	Loose engine mounting screws	Tighten engine mounting screws

\*Centrifuge working hours counter is checked in the settings menu (Hour meter)

## 11. Disinfection and cleaning

Rotor and accessories must be cleaned once a week to avoid corrosion and changes to material.

**Attention! Disconnect the centrifuge from mains power supply, remove rotor, and clean it separately. Before any cleaning or decontamination process other than that is recommended by the manufacturer, the user must ask the manufacturer if the planned process does not damage the device.**

The rotor, rotor chamber, and the outside of the centrifuge should be cleaned with a soaked gauze or soft cloth. For cleaning, water with soap or other water-soluble mild detergent shall be used. To disinfect use an alcohol-based disinfectant (70% isopropanol/water mixture). If corrosive, toxic or radioactive liquids or pathogenic bacteria are spilled in the rotor or its chamber, the centrifuge must be decontaminated thoroughly.

Rotors are high-quality assemblies which withstand extreme stresses. This stability can be impaired by aggressive chemicals. **Aggressive chemicals can cause** corrosion potentially breaking them and damaging the inside parts, leading to electrical failure.

Damage from UV radiation:

Protective coating of the centrifuge protects it from the sun's harmful UV rays.

However, do not use beta, gamma, or any other high-energy radiation for disinfection.

## 12. Transportation and storage

To ensure safety during transportation, the equipment should be packed in the original manufacturers packaging or similar packaging substitute.

Equipment can be transported in any kind of closed transport; make sure that equipment is tightly fixed and transported according to transportation regulations. Equipment should be stored in original manufacturers packaging in dry room with humidity not more than 80% and temperature range of +10 C up to +40 C. It is not recommended to store the equipment for more than 36 month.

## 13. Warranty statements

- The warranty applies to 24-month period from the date of purchasing.
- Malfunctions that occur due to a fault of the manufacturer, during warranty period, are repaired free of charge.
- Warranty is not valid in the following cases:
  - If the serial number label of the manufacturer is damaged.
  - If damage occurs because of
  - incorrect operation, transportation, or storage.
- These documents are necessary if applying for warranty repair:
  - User manual with serial number of the machine.
  - Officially signed report describing the reasons and circumstances of the equipment malfunctions.
- Warranty repair can be performed only if the equipment is delivered in the original manufacturer's packaging or equally safe packaging.
- Before returning for warranty repair ensure that the device is fully cleaned, decontaminated, and does not present any kind of health risk to our staff.
- If the above warranty requirements are not met, repair charges are applied to the customer.
- For all further questions concerning usage and maintenance, please contact the manufacturer or product vendor.

## 14. Disposal

For the disposal of the centrifuge mind the regulations in your country.

## 15. Labelling

### **Centrifuge labelling contains:**

- name of the manufacturer;
- manufacturer's contacts;
- name of the medical device;
- the number of the medical device according to the manufacturer's numbering system (8 characters: XX Y ZZZZ T), where:  
XX - year of production, Y - quarter of production, ZZZZ - batch number, T - centrifuge modification;
- rated mains voltage, Hz frequency; maximum rated power; maximal consumable power;
- the Eurasian Conformity mark;
- CE marking;
- warning symbols;
- RoHS WEEE Recycling Sign.

**Packaging labeling contains:**

- name of the medical device;
- the number of the medical device according to the numbering system of the manufacturer;
- Eurasian Conformity mark;
- Country of Origin;
- handling signs ("Top", "Keep away from moisture");
- name of the manufacturer;
- contacts of the manufacturer.

## 16. Certificate of approval

Centrifuge CM-8S\_\_\_No\_\_\_\_\_ has been inspected for the technical conditions and meets all regulations necessary for this class of device.

Quality control person \_\_\_\_\_  
(name) (signature)

Date of manufacture \_\_\_\_\_

Place for stamp

## 17. Certificate of sale

Organization \_\_\_\_\_

Address \_\_\_\_\_

Phone \_\_\_\_\_

Vendor \_\_\_\_\_  
(name) (signature)

Date of sale \_\_\_\_\_

Manufacturers' address  
Bukultu street 7b, Riga LV-1005  
Tel: (+371) 6755 8743  
(+371) 6755 1934  
E-mail: info@elmi-tech.com  
WWW: www.elmi-tech.com

## ELMI PRODUCT LIST

### CENTRIFUGES BENCHTOP CENTRIFUGES

#### **CM-7S Plus**

3500 RPM desk-top centrifuge for test tubes 10-50 ml.

#### **CM-8S**

4500 RPM desk top centrifuge for operation with centrifuge tubes up to 50 ml.

### FUGAMIX CENTRIFUGE MIXERS

The Fugamix family of products is a unique series, engineered to combine the two most frequently used functions in the laboratory – mixing and spin down – in a single machine.

#### **ELMI CM-50M Fugamix**

Features: 15 000 RPM, 9 levels of mixing and 3 automated spin-mix-spin programs for sample preparation.

#### **ELMI CM-50MP Fugamix Miniprep-Master**

Features: 13500 RPM, Fast Mini Prep, Custom Spin-Mix-Spin Programs

#### **ELMI CM70M-09M Fugamix**

Features: 9000 RPM, 6 levels of mixing, automated spin-mix-spin programs for sample preparation

#### **ELMI CM70M-07M Fugamix**

Features: 7000 RPM, 6 levels of mixing, automated spin-mix-spin programs for sample preparation

### MICROCENTRIFUGES

#### **CM-50**

This mighty mini-centrifuge spins down samples at 15,000 RPM, yet runs on only 24 volts, and fits easily into the smallest spaces.

### DIGITAL THERMO MICROPLATE SHAKERS

**ELMI DTS-4** Digital Thermo Shaker for 4 Micro Plates

**ELMI DTS-2** Digital Thermo Shaker for 2 Micro Plates

### MAGNETIC STIRRER

**ELMI MS-01** 4 Place Magnetic Stirrer

### SKY LINE DIGITAL ORBITAL SHAKERS

**ELMI DOS-20S** Digital Orbital Shaker 20mm Amplitude with Small Platform

**ELMI DOS-20M** Digital Orbital Shaker 20mm Amplitude with Medium Platform

**ELMI DOS-10L** Digital Orbital Shaker 10mm Amplitude with Large Platform

**ELMI DRS-12** Digital Rocking Shaker with 12 Degree Rotation – Rubber Mat Included

**ELMI DOS-20L** Digital Orbital Shaker 20mm Amplitude with Large Platform

**ELMI DOS-10M** Digital Orbital Shaker 10mm Amplitude with Medium Platform

### SKY LINE ANALOG ORBITAL SHAKERS

**ELMI S-3.02 20L** Analog Orbital Shaker 20mm Amplitude with Large Platform

**ELMI S-3.02 20M** Analog Orbital Shaker 20mm Amplitude with Medium Platform

**ELMI S-3.02 10M** Analog Orbital Shaker 10mm Amplitude with Medium Platform

**ELMI S-3.02 10L** Analog Orbital Shaker 10mm Amplitude with Large Platform

### INTELLI MIXERS

**ELMI RM-2L** Intelli-mixer Large - Includes Mix Rack

**ELMI RM-2M** Intelli-mixer Medium - Includes Mix Rack

**ELMI RM-2S** Intelli-mixer Small - Includes Mix Rack

### VORTEX

**ELMI Vortex V-3** with adapter for 0.5 ml to 50 ml volume tubes

### WATER BATH THERMOSTATS

**ELMI TW-2.03** Circulating Water Bath with 8.5 l Plastic Tank, from room temp. to +90°C

**ELMI TW-2.02** Circulating Water Bath with 8.5 l Stainless Steel Tank, from room temp. to +100°C

[www.elmi-tech.com](http://www.elmi-tech.com)

Video instruction is available!

[youtube.com/elmimedia/videos](https://youtube.com/elmimedia/videos)

Video instruction is available!

[youtube.com/elmimedia/videos](https://youtube.com/elmimedia/videos)